



Ionic Oasis®, LLC
Detox Foot Bath

**Best Detox
On The Market**

Eliminates Toxins Naturally



Pg Table of Contents:

1. Ionic Technology Explained
3. Benefits of Detox Foot Baths
4. The Need for Mineral Supplements
5. Disclaimer
6. Frequently Asked Questions
12. About Ionic Oasis Foot Baths

IONIC DETOX TECHNOLOGY EXPLAINED

People love the effects of ionic Detox foot baths; they find themselves more relaxed and relieved of tension and stress. But what is an ionic foot bath and how does it work? One example would be to go to a foot spa session provider and experience a Detox foot bath session for an average fee of about \$30-\$75. Many chiropractors, Physicians, beauty salons some health food stores provide this service. Or you may purchase your own unit.

But what happens during the ionic foot bath session? I will try to explain as simply as possible given that this technology is cutting edge and it would require a lengthy explanation. So I'll stick to the basic facts so as not to confuse the issue.

A Detox Foot Bath operates through a process called electrolysis. This is done by generating the proper amount of DC current in the foot bath water causing the molecules of H₂O to divide producing negative and positive charged particles. These charged particles are called ions. Once the negative ions are present in the water, the body absorbs these ions through a process called osmosis.

So what is osmosis? A simple explanation of osmosis for the application of electrolysis would be the method in which a charged particle can pass through a membrane from a higher concentration to a lower concentration. In other words this is the method in which the ions in your body may become equalized when exposed to an Ionic foot bath. Because people have more positive than negative ions an Ionic foot treatment will theoretically equalize the amount of negative and positive ions in the body. If this is accomplished the person will be in what is called a state of [CHI]. Another term for the word chi when referring to the body would be a state of total harmony and balance. It is an ancient philosophy that when the body is in a state of Chi, this is when the body can best heal itself. A detoxified body is a healthy body and a healthy body can heal itself.

So what is a negative ion? A negative ion is an atom that has lost or gained one or more electrons causing it to become negatively charged. This negative ion is smaller than most molecules which allow it to be easily absorbed by the body through the process of osmosis. Negative ions are produced by the ion field that is present by placing an electrically charged array of electrodes into the water during a foot bath session.

The Ionic Oasis Foot Spa generates millions of these negative ions during the foot bath session. This causes the body to absorb as many negative ions as possible. You could compare this process to a fully charged battery. When it has been fully charged it can not absorb anymore energy because the battery has equalized the negative and positive properties contained in it, and has absorbed as much energy as it can store.

The foreign matter in a person's body such as free radicals, harmful chemicals, heavy metals, parasites, and an endless list of other harmful agents cause many health concerns. It can be very beneficial to ones health to be able to reduce or eliminate these substances. These impurities contain a positive charge, and the negative ions are seeking an opposing or in this instance a positive molecule or particle to attach to. Once attached to these foreign substances they are neutralized then the body can readily identify them and eliminate them naturally through its own

natural processes. Some of the impurities are transferred through the feet during the Detox foot bath session however; once a person has absorbed such a high amount of negative ions their body begins to eliminate these toxins naturally over a period of days.

HOW THE FOOTBATH WORKS

The array goes into the water with the feet, and the power supply delivers a small amount of direct current into the array. When used for feet, this array in a container of water is the ionic foot bath. The electrical charge causes the molecules of the metals within the array, in combination with the water and other elements in the water, to divide thus generating positive and negative ions. This also produces small amounts of oxygen and hydrogen to be released into the air. While you are being treated at the foot spa by the ionic foot bath this is part of the chemistry that takes place when you divide the molecule H₂O.

When experiencing the Detox foot bath, you will see the water change color as the electrical current, water, and magnetic field interact with your feet. While the water would change color by itself, the colors are noticeably different with feet in the foot spa and from one person compared to another person.

BENEFITS AND RISKS OF USING IONIC FOOTBATHS

- **Benefit:** Helps the body **naturally** eliminate harmful toxins and heavy metals and other harmful things compromising good health.
- **Risk:** In the process of natural elimination of toxins it may eliminate essential minerals needed to maintain wellness.
- **Solution:** Talk with your doctor about taking a high quality mineral Supplement while undergoing detoxification.

Note* Before you start any new health regiment check with your Healthcare Provider. They know you and your health condition better then anyone.

We are not medical providers and do not offer medical advice. We are only providing what we consider to be useful information conveyed to us by our customers. This may help you to make an informed decision regarding your personal health prior to purchasing or using of one of our footbaths.

When speaking with your physician here are a few of the questions you may want to ask him/her:

- Do I have an underlying health condition that could cause my body to be over taxed by using a detox footbath?
- Am I on any medications which could be depleting essential minerals in my system?
- Should I take a mineral supplement while using a detox footbath and what should I take?

Ionic Footbaths help the body eliminate harmful heavy metals and toxins from your system, however; it does not discriminate between the good and bad minerals, and may eliminate essential minerals your body needs for good health. Therefore, we highly recommend everyone who is engaged in doing foot bath sessions should speak with their physician about taking a high quality mineral supplement.

TRY THIS...

To prove to yourself or someone else that the footbath is working try using pH test strips prior to your footbath and after. Do this after every session. As your body becomes detoxified your alkaline level will increase.

MINERAL SUPPLEMENT SOLUTION:

Judith Cutler has purchased one of our footbaths for use in her nutritional practice in Southern California and she has developed a special mineral formula supplement made specifically for use with of the Ionic Oasis® foot baths. We are excited about this venture and hope it will be of beneficial use to our customers.

Bio Tech Solutions would be more than happy to do a hair ... to Bio Tech Solutions would be more than happy to assist your needs.

Bio Tech has a [Special Offer](#) for Ionic Oasis® Customers.

SPECIAL OFFER TO IONIC OASIS® FOOT CLEANSE PURCHASERS

You have made an important contribution to your health and well being by committing yourself to a regimen of Ionic Oasis® Foot Cleanses to rid your body of the effects of heavy metals, chemicals and environmental toxins...Well Done!

Here at BioTech Solutions, we use the Ionic Oasis® in our nutritional practice in Southern California and the results are great!

To assist you in your quest for optimum health, we are offering you two services to enhance your progress.

1) An Ionic Mineral Formula – This is a one month formula containing the highest quality ionic minerals that your body needs in order to pull the heavy metals out of your body. It also replaces the minerals that are lost through the detoxification process. It is a ½ oz. dropper bottle of specific ionic minerals. You simply put 20 drops of these minerals in a 6 oz. glass of good water and drink it after each cleanse. It gives you an adequate supply of minerals for 3 cleanses per week (one every other day.) The cost is \$25 per bottle, plus shipping and tax.

2) A Specialized Hair Analysis - This analysis can determine your current state of health and monitor or maximize your detoxification results. We're offering this service at a **Special Discount to Ionic Oasis® Customers!**

Through this test, we can identify the presence of existing heavy metals, chemicals & toxins in your body currently and help to monitor your progress as you continue your cleanse. The cost for this test to Ionic Oasis® customers is \$50 . It is normally a \$75 analysis.

We have been doing specialized hair analyses and customized nutritional formulas since 1989...It's our specialty and we have helped clients all across the country.

If you would like help with your detoxification *or* if you have a specific body problem that you would like to address nutritionally, please contact us at the phone number or e-mail address listed below. We can help!

Judith Cutler
BioTech Solutions
biotheta@yahoo.com
(818) 249-4122

IONIC FOOT BATH DISCLAIMER

The Ionic Foot Bath Detoxification foot spa should not be used by individuals, who are epileptic or hemophiliac, or those with a pacemaker, implanted organs or on blood thinners. It should also not be used by children under age 4, or during pregnancy or lactation. There may be other situations in which you should not use this unit. Please consult your health care provider if you are uncertain.

This product is not evaluated by the Food and Drug Administration. Therefore, these statements contained herein are for informational purposes only.

This product is not intended to diagnose, treat, cure or prevent disease. Ionic Oasis® LLC makes no medical claims or statements of validation. Always consult with your professional Health Care provider before beginning any health care regimen.

The array goes into the water with the hands or feet, and the power supply delivers a small amount of direct current into the array. This causes the metals within the array, in combination with water and salt, to generate positive and negative charged ions by separating oxygen and hydrogen in the water.

WARNING: Never operate any electrical or electronic devices while you are undergoing your foot bath session.

FREQUENTLY ASKED QUESTIONS

[Can anyone use the foot bath?](#)

[Is this an approved FDA product?](#)

[Is there a warranty?](#)

[Is this a medical device, can it cure or treat my illness?](#)

[Why does the water change color without feet in the water?](#)

[How often can I do my sessions?](#)

[Why can't I do daily sessions?](#)

[How does my body detoxify in the water?](#)

[Why does the water color change?](#)

[Why does the water color vary from treatment to treatment?](#)

[Can I touch the array during treatment?](#)

[Can I use the FootBath while taking prescription drugs?](#)

[Should I remove my jewelry?](#)

[Is it safe people who have metal plates in their body?](#)

[Why do some people feel tired after a treatment?](#)

[What is the function of the array?](#)

[What are the bubbles coming from the array?](#)

[How do I preserve and clean the array?](#)

[How long will the array last?](#)

[How can I tell when to replace the arrays?](#)

[Why do my feet tingle?](#)

[How much are the array's rods.](#)

[What brand is this bath?](#)

[Does this create negative and positive ions?](#)

[Do your units alternate polarity like most of the higher cost units?](#)

[What happens if the array element tips over in the water?](#)

[Do you make these units yourself?](#)

[Does the unit have an automatic timer to shut off like other units?](#)

[What is the difference between the 3/8 inch to the 1/2 inch?](#)

[How does the array keep from falling over in the water?](#)

[What makes your unit better than some of the others out there?](#)

[Does this unit offer the same level and quality of Detoxification?](#)

[Are your units dangerous and at risk of overheating?](#)

[Do you personally manufacture these units?](#)

[Do you offer wholesale pricing?](#)

[What makes the PRO model better?](#)

[Can the PRO unit stand up to multiple uses?](#)

[Can the PRO unit stand up to several years of moderate to heavy use?](#)

CAN ANYONE USE THE FOOT BATH?

- The Ionic Foot Bath Detoxify foot spa should not be used by individuals, who are epileptic or hemophiliac, or those with a pacemaker, implanted organs or on blood thinners. It should also not be used by children under age 4, or during pregnancy or lactation. There may be other situations in which you should not use this unit. Please consult your health care provider if you are uncertain.

IS THIS AN APPROVED FDA PRODUCT?

- Detox Foot Baths have not, nor has the technology of Detox Foot Baths been evaluated or approved by the (FDA) Food and Drug Administration. Therefore, these statements contained herein are for informational purposes only.

IS THERE A WARRANTY?

- Warranty of Array: Our arrays carry a five (5) year warranty to the original owner against defects in materials and workmanship. During the first six (6) months, replacement parts and shipping are free. After six (6) months you pay only for shipping. A receipt or proof of payment is required. Warranty does not cover normal wear of replaceable elements. Damage due to abuse is not covered.
- Warranty of Power Supply: Our power units carry a two (2) year manufacturer's warranty.

IS THIS A MEDICAL DEVICE, CAN IT CURE OR TREAT MY ILLNESS?

- This product **is not intended** to diagnose, treat, cure or prevent disease. Ionic Oasis® makes no medical claims and makes no statements of validation. Always consult with your professional Health Care provider before beginning any health care regiment.

WHY DOES THE WATER CHANGE COLOR WITHOUT FEET IN THE WATER?

- If you operate the foot bath without your feet in the water, it will turn the water different colors, yellow to orange to brown because the element usage, minerals, and other impurities in the tap water will precipitate out. However; you can see the difference between this and a footbath session. **Please note** that if the tap water is particularly impure, the difference in water color between feet and no feet will be less noticeable, just like adding a handful of dirt to already dirty water. The fact that the water changes color is evidence that the array is working properly.

HOW OFTEN CAN I DO MY SESSIONS

- Between the ages of 10 and 65, every other day not to exceed 3 times per week.
- Under 10 and over 65 the recommendation is no more then twice a week.

WHY CAN'T I DO DAILY SESSIONS?

- Daily treatments may create a healing crisis as the body naturally eliminates toxins. It may make you feel excessive fatigue.

HOW DOES MY BODY DETOXYFY IN THE WATER?

- An energized and balanced body will naturally be able to detoxify at a higher rate.

WHY DOES THE WATER COLOR CHANGE?

- It is due to the interaction of the metals from the array, sea salt, and minerals and impurities contained in the water along with the individuals contribution, creating an ionic field with the water.
- The water will change color regardless of whether you have an individuals feet in the water or not, because the metallurgic reaction of the array and the polarization of the solids and minerals in the water. However; there will normally be 2-3 times the cellular debris within the water when the session is combined with a individual.

WHY DOES THE WATER COLOR VARY FROM TREATMENT TO TREATMENT?

- Variations occur due to each individual's energy state and toxicity levels as toxins are re-ionized at the molecular level and react with the metals of the array.

CAN I TOUCH THE ARRAY DURING TREATMENT?

- Absolutely, this will not alter the session in any way.
- **WARNING:** Never operate any other electrical or electronic devices while you are undergoing your foot bath sessions.

CAN I USE THE FOOTBATH WHILE TAKING PRESCRIPTION DRUGS?

- Most doctor prescribed medications will not be affected by use of the footbath, but the decision is best left up to you and your physician. We recommend consulting your physician prior to use.

SHOULD I REMOVE MY JEWELRY?

- Any Jewelry that will be in the water should be removed.

IS IT SAFE PEOPLE WHO HAVE METAL PLATES IN THEIR BODY?

- It is best to speak with your health care provider regarding your concerns with metal implants. The other concerns are noted in our [disclaimer](#).

WHY DO SOME PEOPLE FEEL TIRED AFTER A TREATMENT?

- Some people may experience fatigue as the body's detoxification process accelerates itself. This will balance itself as the body expels the toxins naturally and regenerates itself.

WHAT IS THE FUNCTION OF THE ARRAY?

- The array is the conductive transfer unit of the direct current into the water which generates the ionic field.

WHAT ARE THE BUBBLES COMING FROM THE ARRAY?

- The bubbles are the separation of the hydrogen and oxygen in the water as the ionic field is being created.

HOW DO I PRESERVE AND CLEAN THE ARRAY?

- Initially you can rinse your array with very warm water and pat dry with a towel. However; as the elements become pitted from use, it may be put into the dishwasher or follow the next step below.
- Dip the array in a solution mixture of 1 part vinegar to 2 quarts warm water for a couple of minutes. Then brush the array to remove access particles and rinse thoroughly with warm water and pat dry.
- If you are using the footbath in a professional setting we would recommend the use of a good disinfectant cleaner rather than the steps noted above.

HOW LONG WILL THE ARRAY LAST?

- The array itself if properly cared for will last indefinitely however; the generating rods in the array will eventually lose their effectiveness after approximately 80-120 sessions for the Model I, MiniOne and 120-160 sessions for the Models II, Model Pro, and MiniTwo and MiniPro depending on the water conditions. They will eventually lose their conductive properties and will need to be replaced.
- You can purchased a new set of 304 stainless steel rods from us at <http://www.ionicoasis.com/ionic-footbath-purchase/>.

HOW CAN I TELL WHEN TO REPLACE THE ARRAYS?

- The 304 stainless pods will be noticeably deteriorated after numerous sessions. You will be able to tell when the water does not change color as it should. NOTE* You can reverse the terminals on the array when it begins to deteriorate. This will make it last another 20-40 sessions due to the fact that the positive electrodes wear out a lot sooner then the negative terminals. Reversing the connections will not affect the operation of the array but only allow the elements to wear more evenly thus making it last longer.

WHY DO MY FEET TINGLE?

- It is normal to feel a slight tingling sensation during a session, but if you experience uncomfortable tingling or itching this is an indication there is too much salt in the water. Try using less salt on your next session. Too much salt will not alter the effects of the session, however; increasing the salt also increases the amperage of the current. Too much salt can result in overheating the power unit causing it to shut down. The power supply will not resume operation until it has fully cooled down usually in about 25 minutes. As per the operating instructions you should begin with 4 1/2 to 5 inches of water and start with 1 scoop of salt. If you prefer you can increase the amount of salt but not to exceed 2 scoops. After you have operated the system several times you will know how to set it up for each use.
- We recommend using permanent marker to mark the water level at 4 1/2 to 5 inches if you are using a tub that is not already marked. This will ensure the same results with each use with the same amount of water and salt added. Check our [Instruction Manual](#) for updated chart information regarding tub sizes, amount of salt to use, and water conditions.

HOW MUCH ARE THE ARRAY'S RODS.

- You can find our rod replacements at <http://www.ionicoasis.com/ionic-footbath-purchase/>

WHAT BRAND IS THIS BATH?

- This is an Ionic Oasis® models manufactured by Ionic Oasis® LLC in Lincoln and Beatrice, NE.

DOES THIS CREATE NEGATIVE AND POSITIVE IONS?

- Yes it does, due to the chemical composition of water. It creates equal amounts of negative and positive ions. It is not possible to create a negative ion without also creating a positive ion simultaneously. Is is the negative ions which are beneficial.

DO YOUR UNITS ALTERNATE POLARITY LIKE MOST OF THE HIGHER COST UNITS?

- No our units do not alternate from positive to negative. All our units produce positive and negative ions in equal quantities, at the same time. It is impossible to generate negative without also generating positive ions simultaneously. We do not know how it would be possible to produce just negative or just positive ions. Maybe one of these other companies could explain this claim. We believe that we should only state the known facts about our product, because we don't want our customers to be misled.

WHAT HAPPENS IF THE ARRAY ELEMENT TIPS OVER IN THE WATER?

- Our units come with a GFCI Adapter which is a circuit breaker which breaks the electric current. The power unit plugs into the Ground Fault Plug Interrupter and the GFCI plugs into the wall. If the array tips over in the water simply pick it up and set it back upright. This would not be a normal incident as our design of the lexan legs should prevent this from happening under normal use.

DO YOU MAKE THESE UNITS YOURSELF? ALSO, IS IT EASY TO SHOCK ONESELF ON THEM?

- Yes, we do manufacture these products.
- No, it is not easy to shock oneself on them. You can in fact touch the array with your hands or feet while it is in operation and not feel any more than a slight tingle at most.
- If the product is used as directed, there is no danger. Not only is the power unit fused, but the GFCI adapter will also insure your safety.

DOES THE UNIT HAVE AN AUTOMATIC TIMER TO SHUT OFF LIKE OTHER UNITS?

- No our units do not have automatic timers.

WHAT IS THE DIFFERENCE BETWEEN THE 3/8 INCH TO THE 1/2 INCH?

- All our units perform very well however the 1/2 inch units will give about 40% more sessions before the elements need to be replaced. Also our model pro and Mini Models will yield the same results in 30 minutes that the model I and model II will yield in 40 minutes.

HOW DOES THE ARRAY KEEP FROM FALLING OVER IN THE WATER?

- The legs on the array are triangular in shape which stabilizes the legs. The legs are made of sturdy lexan plastic, with 2 the cross bars on top. Between the shape of the legs, and materials used the array, and the weight of the array, is very sturdy practically indestructible. When the array is in the foot bath it does not have a tendency to tip. If by chance it would get tipped over in the water it will not hurt anything you simply pick it up and set it back up.

WHAT MAKES YOUR UNIT BETTER THAN SOME OF THE OTHERS OUT THERE AND WHY DO YOU FEEL OTHERS CHARGE SO MUCH FOR THEIR DIGITAL AND BETTER VISUALLY PACKAGED UNITS?

- I can not say for sure why other company's charge higher prices for their products. I would like to think it is to re-coop their overhead: marketing, name brand, research, etc. We do not have this type of overhead which is why we can charge much less for our products.

DOES THIS UNIT OFFER THE SAME LEVEL AND QUALITY OF DETOXIFICATION AS OTHER MORE HIGH PRICED UNITS OUT THERE SAY THAT THEY DO?

- Yes it does. A Chiropractor in IN called to commend us on our product as he was paying \$45.00 for his Detox sessions at a clinic, and while he was skeptical if ours would do as promised, he said ours did as well if not better than the high priced model he used at a clinic and thanked us for making an affordable product. This was not an email but a personal call. I think his call speaks volumes, as he didn't have to call or email to tell us, but he did.

ARE YOUR UNITS DANGEROUS AND AT RISK OF OVERHEATING?

- No they are not dangerous. The power supply is Thermal Overload Protected and comes with a 2 year manufacturer warranty. The only way to overload the power supply is by using too much salt. When this happens the power supply will shutdown for about 25 minutes. Once it has cooled you are ready to begin your next session.

DO YOU PERSONALLY MANUFACTURE THESE UNITS?

- We manufacture the array's ourselves. We purchase the power supplies, pans, and the salt, from other distributors.

DO YOU OFFER WHOLESALE PRICING IF I WANTED TO BUY IN A QUANTITY OF 10-20 PCS.?

- The volume wouldn't be a problem, we do have a Dealership Agreement for anyone interested in selling our products, please contact OasisTeam@ionicoasis.com

WHAT MAKES THE PRO MODEL BETTER?

- The PRO model is geared more toward users who wish to use their footbaths in salons, clinics, in-home spa sessions etc. With the Pro Model you achieve the same results with 30 minutes sessions and the ½" elements last 40% longer.

CAN THE PRO UNIT STAND UP TO MULTIPLE USES?

- Our models have been tested non stop for a period of 48 hours at full power. So multi use will not be a problem.

CAN THE PRO UNIT STAND UP TO SEVERAL YEARS OF MODERATE TO HEAVY USE?

- Yes, as a matter of fact we offer a 5 year warranty on the array. Actually, all of our models carry a 5 year warranty. However; warranty does not cover the elements. It is up to the purchaser to replace elements after 80-120 or 120-160 sessions depending on the model you purchase.

ABOUT IONIC FOOT BATHS:

This system is designed to be simple reliable safe to use, easy to maintain, and built to last. In fact this system out performs most other units on the market at a much more affordable price.

Our 3/8" elements will last 80 to 120 and the 1/2" elements will last 120 to 160 sessions before the generating elements need to be replaced.

The Ionic Foot Baths Detoxification foot spa should not be used by individuals, who are epileptic or hemophiliac, or those with a pacemaker, implanted organs or on blood thinners. It should also not be used by children under age 4, or during pregnancy or lactation. There maybe other situations in which you should not use this unit. Please consult your health care provider before using our Foot Spa if you are uncertain.

This product has not evaluated by the Food and Drug Administration. Therefore, these statements contained herein are for informational purposes only.

This product is not intended to diagnose, treat, cure or prevent disease. Ionic Oasis® makes no medical claims and makes no statements of validation. Always consult with your professional Health Care provider before beginning any health regimen.

ABOUT IONIC OASIS® DETOX FOOT SPA ARRAYS:

- All of our units are made with the highest quality materials for a lifetime of trouble-free use.
- The screw threads are cold formed in the aluminum, not tapped, as is common. This results in a much stronger and more durable screw thread.
- The holders for the electrodes are made of 60-61 Aircraft Quality Aluminum.
- The elements are 18-8, 304 Stainless Steel. This is an Austenitic, Chrome-nickel, food grade stainless used in the food industry and many other heavy industrial applications.
- The legs of the array are made of lexan plastic. This is a very durable material used in the aerospace industry for items such as aircraft canopies, windscreens and other windows. It is also one of the most common products used to make bullet-resistant glass.
- And the portable GFCI adapter will protect you from any unforeseen risk of electrical shock. Not all GFCI's are not created equal. Our high quality GFCI detects changes in current as low as 4 MilliAmps allowing it to trip (turn off) instantly under the most minute change in current. Assuring the user a high level of protection and peace of mind that our unit is safe under all conditions.

WARRANTY OF ARRAY:

Our arrays carry a five (5) year warranty to the original owner on defects in materials and workmanship. A receipt or proof of payment is required. Warranty does not cover normal wear of replaceable elements. Damage due to abuse is not covered.

WARRANTY OF POWER SUPPLY:

Our power units carry a two (2) year manufacturer's warranty.

NOTE* Operating the foot bath without feet in the water, will turn the water different colors; from yellow to orange to brown because the minerals and other impurities in the tap water will precipitate out. However; you can see the difference between this and a footbath session. Please note that if the tap water is particularly impure, the difference in water color between feet and no feet will be less noticeable, just like adding a handful of dirt to already dirty water. The fact that the water changes color is evidence that the array is working properly.